


## **DENSITY CONTROL KIT** complementary Anchor strand 1 and 2

Dynamic Biosensors GmbH & Inc.  
DC-0 v5.1



## Key Features

- Complementary **Anchor strand 1** and **complementary Anchor strand 2** for low density generation of **helix<sup>®</sup> Adapter Chip Spot 1** and **Spot 2**.
- Compatible with **helix<sup>®</sup> Adapter Chip**.
- Ideal for **MIX&RUN** sample preparation.
- This kit enables to modulate the ligand density on the sensor surface.

## helix<sup>®</sup> Adapter Chip Overview

2 spots with 2 different anchor sequences for DNA-encoded addressing.



## Product Description

Order Number: **DC-0**

Table 1. Contents and Storage Information

Material	Cap	Concentration	Amount	Buffer	Storage
<b>cAnchor strand 1</b>	Red	200 nM	5 x 200 $\mu$ L	TE40 <sup>[1]</sup>	-20°C
<b>cAnchor strand 2</b>	Red	200 nM	5 x 200 $\mu$ L	TE40 <sup>[1]</sup>	-20°C

For research use only.

This product has a limited shelf life, please see expiry date on label.

To avoid many freeze thaw cycles please aliquot the nanolever.


## Preparation | MIX&RUN

In-solution hybridization of adapter and ligand strands:

1. Mix **Adapter strand 1 - Ra** (400 nM) and conjugated **Ligand strand** (500 nM) at 1:1 ratio (v/v).
2. Incubate the solution of step 1 at **RT** at **600 rpm** for **30 min** to ensure complete hybridization.
3. Mix solution of step 2 and **Adapter strand 2 - Ra - lfs** (200 nM) at 1:1 ratio (v/v).

Solution is ready to use for biochip functionalization.

Stability of the solution is related to the stability of the ligand molecule.




## Preparation | Ligand Density Variation

Easy modulation of ligand density on the sensor surface:

4. Mix **cAnchor strand 1** and **cAnchor strand 2** at 1:1 ratio (v/v).
5. Mix solution of step 3 and solution of step 4 at the ratio of desired on-chip density, e.g.:
  - a. 1:1 ratio (v/v) for a 50% ligand density
  - b. 1:2 ratio (v/v) for a 30% ligand density
  - c. 1:9 ratio (v/v) for a 10% ligand density

Solution is ready to use for biochip functionalization.



## Example

Required volume for 3 functionalizations: **100 µL** with a final **ligand density** of **10 %**.

Vial 3	Vial 4	
<b>Adapter strand 1 - Ra</b> with ligand (100 nM) <b>Adapter strand 2 - Ra - lfs</b> (100 nM)	<b>cAnchor strand 1</b> (200 nM)	<b>cAnchor strand 2</b> (200 nM)
10 µL	45 µL	45 µL

Mix vial 3 and vial 4 to obtain 100 µL of ready-to-use DNA solution.

## Useful Order Numbers

Table 2. Order Numbers

Product Name	Comment	Order No
<b>helix® Adapter Chip</b>	Chip with 2 detection spots	ADP-48-2-0
<b>helix® Amine Coupling Kit 1</b>	For five individual conjugation reactions	HK-NHS-1
<b>Adapter strand 1 - Ra</b>	400 nM	AS-1-Ra
<b>Adapter strand 2 - Ra - lfs</b>	200 nM, prehybridized with Ligand-free strand	AS-2-Ra-lfs

## Contact

**Dynamic Biosensors GmbH**  
Perchtinger Str. 8/10  
81379 Munich  
Germany

**Dynamic Biosensors, Inc.**  
300 Trade Center, Suite 1400  
Woburn, MA 01801  
USA

**Order Information** [order@dynamic-biosensors.com](mailto:order@dynamic-biosensors.com)  
**Technical Support** [support@dynamic-biosensors.com](mailto:support@dynamic-biosensors.com)

[www.dynamic-biosensors.com](http://www.dynamic-biosensors.com)

Instruments and chips are engineered and manufactured in Germany.  
©2024 Dynamic Biosensors GmbH | Dynamic Biosensors, Inc. All rights reserved.