

Biophysical Analysis of Cas9-DNA Interactions and Enzymatic Activity



with the switchSENSE® Biosensor Platform

Felix Kröner, Prayanka Rajendran and Ulrich Rant; Dynamic Biosensors GmbH, Martinsried, Germany

The switchSENSE® Principle: Electro-Switchable DNA Nanolevers

DNA nanolevers are electrically actuated at high-frequency on microelectrodes, while their orientation is monitored by time-resolved single photon counting. The binding of analyte molecules slows the switching dynamics in a characteristic way.


Nature Commun. 4:2099 (2013) | Bioanal. Rev. 4 (2) 97-114 (2012) | JACS 134, 15225 (2012) | PNAS 107, 1397 (2010) | JACS 132, 7935 (2010)



DRX & DRX²
Limit of detection 10 fM
Dissociation constant 50 fM – 1 mM
Association rate constant 1E3 – 1E8 /Ms
Dissociation rate constant 1E-6 – 1E0 /s
Hydrodynamic diameter accuracy of 0.1 nm
Temperature 8° – 75°C




switchSENSE® biochip
4 flow channels
20 detection spots á 0.01 mm²
Multiplexing via DNA encoded addressing




High Sensitivity | Kinetics and Affinity | Size and Conformation | Cooperativity and Avidity | Thermodynamics

Investigating Cas9 with the switchSENSE® platform



dCas9 vs. Cas9 vs. apoCas9 kinetics




↑ Binding curves of dCas9 and Cas9 (c = 20 nM) show weak off-target binding. Their binding kinetics are monitored through fluorescent proximity sensing. Cas9 shows stronger fluorescence decrease than dCas9, indicating DNA cleavage due to fluorophore removal. We find an on-rate for dCas9 and Cas9 of $k_{on} \approx 1 \cdot 10^6 \text{ M}^{-1}\text{s}^{-1}$.

↑ apoCas9 (c = 20 nM) binds DNA sequence independent as found by Doudna et al. (2014)¹. We find PAM independent of apoCas9 (data not shown), with an on-rate, $k_{on} = 2.11 \pm 0.01 \cdot 10^5 \text{ M}^{-1}\text{s}^{-1}$. We hypothesize unspecific binding of the free gRNA binding site of apoCas9 to NAs.

¹Sternberg, S. H., Redding, S., Jin, M., Greene, E. C., & Doudna, J. A. (2014). Nature, 507(7490), 62–67. <http://doi.org/10.1038/nature13011>


So far unidentified fast on-rates with switchSENSE®



← Mass-transport limitation: Conventional biosensors with high-ligand densities and non-optimal flow conditions are prone to measurement artifacts like mass-transport limitations and rebinding. switchSENSE® reveals faster on-rates than previously identified for dCas9/Cas9. High flow rates reveal faster on-rates indicating mass transport limitation. The found on-rate, $k_{on} = 1.5 \cdot 10^6 \text{ M}^{-1}\text{s}^{-1}$ is about 50 times faster than the on-rate found by Richardson et al. through BLI measurements.


²Richardson et al. Nat Biotech 34 (3) 339 (2016)

Exchangeable target sequence



Positioning the target sequence above the fluorophore provides an easy and fast way for exchanging the target sequence. The binding of Cas9 and dCas9 is very well observable and the found on-rates agree with the above presented setup. The use of molecular dynamics measurements allows detection of the cleavage of DNA by Cas9. This presents an efficient setup for investigating sequence dependent interactions of any kind of nucleic acids with Cas9.

dCas9 vs. Cas9 cleavage activity



↑ Cleavage efficiency and specificity of dCas9 and Cas9. Comparing fluorescence levels before and after removal of dCas9/Cas9 reveals only cleavage activity of Cas9. Cas9 drastically reduces the fluorescence of the target strands which directly correlates to DNA cleavage. We also find weak off-target cleavage by Cas9 ($\approx 5\%$ of the total cleavage activity).