

10x Buffer TE40 pH 7.4

After dilution this buffer is intended for use as a running buffer for the **heliX**[®] and DRX series.

Product Description

Order Number **BU-TE-40-10**

Material	Amount	Storage	Comment
10x Buffer TE40 pH 7.4 (100 mM Tris, 400 mM NaCl, 500 µM EDTA, 500 µM EGTA and 0.5 % Tween20; 0.2 µm sterile filtered)	50 mL	2-8°C	

For *in vitro* use only.

This product has a limited stability, please see expiry date on label.

Preparation

Dilute 10x Buffer TE40 pH 7.4 with deionized, 0.2 µm filtered water by mixing the complete solution (50 mL) with 450 mL ultrapure water. After dilution TE40 Buffer is ready for use (10 mM Tris, 40 mM NaCl, 50 µM EDTA, 50 µM EGTA and 0.05 % Tween20).

The diluted buffer should be stored at 2-8°C.

Contact

Dynamic Biosensors GmbH

Perchtinger Str. 8/10
81379 Munich
Germany

Phone: +49 89 89 74 544 0

Dynamic Biosensors Inc.

300 Trade Center, Suite 1400
Woburn, MA 01801
USA

Phone: +1 781 404 6126

Order Information


Phone: +49 89 89 74 544 0

Email: order@dynamic-biosensors.com

Technical Support

Phone: +49 89 89 74 544 66

Email: support@dynamic-biosensors.com



switchSENSE® is a proprietary measurement technology by Dynamic Biosensors GmbH. Instruments and biochips are engineered and manufactured in Germany.

©2023 Dynamic Biosensors GmbH | Dynamic Biosensors Inc. All rights reserved.