

proFIRE® Antibody Oligo Conjugation Kit

Functionalization of DNA via amines (-NH₂)




Key Features

- Coupling of antibodies with primary amines (e.g. NH₂-terminus, lysines) to modified DNA in a reaction tube
- Sufficient reagents for 3 reactions plus one control oligo
- Convenient standard chemistry
- Applicable for any antibody formats
- Coupling of multiple antibodies can be performed simultaneously
- Yields > 95 % pure antibody-DNA conjugate with controlled quality of your product
- With any DNA sequence and length up to 150 bases feasible

Workflow Overview

3-Step Conjugation Workflow (in-vitro)



1. DNA Modification

The DNA is activated with amine reactive groups.

2. Antibody Conjugation

After incubation the excess linker is removed by a spin column. The antibody is added to the functionalized DNA and incubated for at least 1 h.

3. Purification

The antibody-DNA conjugate is purified using the **proFIRE®** system.

4. Ready-to-use fractions

The fractions with antibody-DNA conjugate are ready for further processing.

Time line: Hands on time < 1 h | Incubation ~ 2 h | Total ~ 3 h

Product Description

Order Number **PF-AB-1**

TABLE 1 | Contents and Storage Information

| Material | Cap | Amount | Storage | Comment |
|--|-------------|------------|---------|---------|
| Antibody Buffer | transparent | 5 x 1.8 mL | -20°C | |
| ddH ₂ O | transparent | 1.5 mL | -20°C | |
| Crosslinker | brown | 4 x | -20°C | |
| Control oligo, 48 bases with 5'-DBCO | blue | 1 x | -20°C | |
| Purification spin column | red | 8 x | 2-8°C | |
| 2.0 mL Reaction tubes for Purification spin column | | 8 x | r.t. | |
| Centrifugal filter unit (30 kDa MWCO) ¹ | | 4 x | r.t. | |
| Centrifugation collection tube | | 8 x | r.t. | |

For *in vitro* use only.

Please check date of expiry on the kit. Products are shipped at ambient temperature.

The kit contains reagents sufficient for 3 conjugations of approx. 50-500 µg antibody each and for one positive control (control oligo). The resin slurry of the Purification spin column contains 0.02 % sodium azide.

¹ See page 10 for order no.

Additional Materials Required

TABLE 2 | Additional Materials.

| Material | Comment |
|-------------------------------------|---|
| Oligo with DBCO modification | We recommend using 3 - 4 nmol DNA (modified with DBCO, HPLC grade) for one reaction |
| Benchtop microcentrifuge | Capable between 1,000 x g and 13,000 x g |
| Vortexer | |
| 1.5 mL reaction tubes | |
| UV-Vis spectroscopy (e.g. Nanodrop) | Concentration determination of the conjugate |

All necessary solutions and buffers are included in the kit.

Important Notes

- Do not use any buffer containing primary amines (i.e. TRIS, glycine) during conjugation process.
- Dithiothreitol (DTT) can be used up to 1 mM during the conjugation process. Do not use 2-Mercaptoethanol or other thiol-based reducing agents during conjugation process. If a reducing agent is necessary, TCEP is recommended up to 1 mM.
- Avoid using partially purified protein samples or protein samples containing carriers (e.g. BSA).
- To get highest reaction yields, the antibody should be dissolved in Antibody Buffer. Buffer exchange is recommended prior to conjugation process¹.
- Before you begin, briefly centrifuge all tubes with blue, brown and transparent caps to ensure that all material is at the bottom of the tubes.

3-Step Conjugation of a Antibody to a Nanolever in a Reaction Tube

Please read the entire protocol before starting and **perform conjugation without interruption**.

TIP: The protocol can be performed simultaneously for multiple coupling reactions.

Note: Each conjugation kit is supplied with a control oligo (a 48 base oligo with a 5'-DBCO). This oligo is included as positive control in order to give the option of confirming the conjugation chemistry is working optimally. Please proceed with the control oligo as with your oligo of choice.

Before you begin: Allow the crosslinker to reach room temperature before use.

I Nanolever Modification

1. Dissolve your oligo of choice in **40 µL Antibody Buffer** prior to use and vortex until solids are completely dissolved and spin down shortly.
Note: Please proceed with the control oligo as with your oligo of choice.
2. Dissolve the crosslinker (brown cap) by adding **100 µL** ddH₂O and vortex until solids are completely dissolved and spin down shortly. **IMPORTANT:** Always use fresh compounds.
3. Add **10 µL** of the freshly prepared linker solution to one DNA aliquot. Discard the remaining linker solution from step 2.
4. Vortex the reactants for 10 sec, spin down and incubate them for **20 minutes** at room temperature. **IMPORTANT:** Do not exceed incubation time as the reaction yield will decrease.
5. In the meantime, equilibrate **two** purification spin columns (red cap) for one coupling reaction:
 - a. Remove column's bottom closure and loosen cap (do not remove cap).
 - b. Place column in a 2.0 mL reaction tube.
 - c. Centrifuge at 1,500 × g for 1 minute to remove the storage solution.
 - d. Add **400 µL of Antibody Buffer** on top of column's resin bed. Centrifuge at 1,500 × g for 1 minute to remove buffer.
 - e. Repeat step d once, discard buffer from the reaction tube. The Purification spin column should be in a dry state now.
6. Sample loading
 - a. Place columns from step 5 in new 1.5 mL reaction tubes.
 - b. Remove cap of spin column number 1 and apply the sample from step 4 to the top of the resin bed.
 - c. Centrifuge at 1,500 x g for 2 min to collect the sample (flow-through).
Discard Purification spin column after use.
 - d. Remove cap of spin column number 2 and apply the sample from step c on top of the resin bed.
 - e. Centrifuge at 1,500 x g for 2 min to collect the sample (flow-through).
Discard Purification spin column after use.

II Antibody Conjugation

7. Add approx. **100 µg** (up to 500 µg) of the antibody (concentration approx. 0.5 – 50 mg/mL) to the sample from step 6. For optimal conditions use a volume of approx. 50 µL.

EXAMPLE: Adjust antibody concentration to 5 mg/mL and use 20 µL for conjugation.

IMPORTANT: Be sure that the storage buffer of the antibody does not contain any primary amines, e.g. TRIS buffers, glycine (please see page 4, Important Notes).

8. Mix the reaction by pipetting up and down and let it react at room temperature for **at least 1 hour**.


IMPORTANT: Do not vortex. If necessary, the reaction can be carried out at 4 °C with a longer reaction time (e.g. overnight).

III proFIRE® Purification


Please refer to the **proFIRE®** User Manual.

9. Perform a purification using the **proFIRE®**. Please make sure that the sample volume is 160 µL.
- If the volume is less than 160 µL, add Antibody Buffer.
 - If it exceeds 160 µL, please perform two subsequent runs.
10. Use the Data Viewer software of the **proFIRE®** to identify which fractions contains pure conjugate.

Example **proFIRE®** chromatogram of a **100 µg** Antibody conjugation to a 48 mer oligo:



Example **proFIRE**® chromatogram of a **500 µg** Antibody conjugation to a 48 mer oligo:




Used buffers: **proFIRE**® Buffer A; **proFIRE**® Buffer B.
 Column: **proFIRE**® column. Flow: 1 mL/min.
 Used program: DNA length 48. Type: 1.

Note:

The crosslinker will be linked to the amino groups (NH₂-groups) of antibodies. Among the four amino acids with side-chain amino groups (glutamine, lysine, arginine, asparagine), only lysine can react. In addition, the N-terminal amino group can be functionalized with the oligo, too. Thus, there are many sites that can be linked via the crosslinker as shown in the diagram below.

Amino acids with reactive NH₂-groups:

- (i) Lysine (K)
- (ii) NH₂-terminus of the antibody



11. Take the recommended fractions out of the fraction collector.
12.
 - a. Option 1: Store fractions between 8 °C and -86 °C as desired.
 - b. Option 2: Proceed with Buffer Exchange and Concentration (see section IV).


IV Optional: Buffer Exchange and Concentration

13.
 - a. Add **500 µL** of the first fraction containing the antibody-DNA conjugate from the **proFIRE®** to the centrifugal filter unit.
Centrifuge at 13,000 x g (up to 14,000 x g) for **5 minutes** and discard flow-through.
 - b. Add the remaining fractions in the same filter unit and repeat the centrifugation step in order to collect all samples in one tube (Please check on page 9: Additional information for the right use of centrifugal filter unit).
 - c. Add **350 µL** of the buffer of choice for buffer exchange and centrifuge at 13,000 x g for **5 minutes**. Discard the flow-through again.
 - d. Add **350 µL** of the buffer of choice for buffer exchange and centrifuge at 13,000 x g for **10 minutes**. Discard the flow-through again.
 - e. To recover the antibody-DNA conjugate, place the centrifugal filter unit upside down in a **new** centrifugal collection tube (provided in the kit).
Spin for **2 minutes** at 1,000 x g to transfer the sample to the tube.
14. Check antibody-DNA conjugate concentration after buffer exchange by using absorbance at 260 nm and the following equation:
$$c (\text{antibody-DNA conjugate}) = A_{260 \text{ nm}} / (\epsilon * d)$$


ϵ = Extinction Coefficient of the DNA
 d = optical path length
(usually $d = 1 \text{ cm}$, please check photometer manual for further information).
15. Store between 8 °C and -86 °C as desired.

Additional Information

I Buffer Exchange and Concentration with Centrifugal Filter Units



1. Take one centrifugal filter unit, add the appropriate volume of buffer in the filter device, and cap it.
2. Place capped filter device into the centrifuge rotor, aligning the cap strap toward the center of the rotor; counterbalance with a similar device.
3. Spin the device at 13,000 x g (or 14,000 x g) for the given time.
4. Remove the flowthrough and repeat the steps 1-3.
5. Remove the assembled device from the centrifuge and separate the filter device from the microcentrifuge tube.
6. To recover the conjugate, place the filter device upside down in a clean centrifugal tube, aligning open cap towards the center of the rotor; counterbalance with a similar device. Spin for 2 minutes at 1,000 x g to transfer the sample from the device to the tube.



Useful Order Numbers

TABLE 3 | Order Numbers.

| Product name | Order Number |
|--|--------------|
| proFIRE® Amine Coupling Kit 1 for proteins (>5 kDa); sufficient for 5 conjugation series | PF-NH2-1 |
| proFIRE® Thiol Coupling Kit 1 for proteins (>5 kDa); sufficient for 5 conjugation series | PF-SH-1 |
| Centrifugal filter unit (3 kDa MWCO), 5 pcs. | CF-003-5 |
| Centrifugal filter unit (10 kDa MWCO), 5 pcs. | CF-010-5 |
| proFIRE® column | PF-CC-1 |
| 10x proFIRE® Buffer A (50 mL) | PF-BU-A-10 |
| 5x proFIRE® Buffer B (50 mL) | PF-BU-B-5 |
| 1x Conjugation Buffer (12 mL) | PF-BU-C-1 |

My Notes

Contact

Dynamic Biosensors GmbH

Perchtinger Str. 8/10
81379 Munich
Germany

Dynamic Biosensors, Inc.

300 Trade Center, Suite 1400
Woburn, MA 01801
USA

Order Information order@dynamic-biosensors.com

Technical Support support@dynamic-biosensors.com

www.dynamic-biosensors.com



Get it on [Google Play](#).
Download on the [App Store](#).

switchSENSE® and **proFIRE®** is a proprietary measurement technology by Dynamic Biosensors GmbH. Instruments and biochips are engineered and manufactured in Germany.

©2023 Dynamic Biosensors GmbH | Dynamic Biosensors, Inc. All rights reserved.