


## 1X BUFFER M PH 6.5

coupling buffer for nanolever conjugation

Dynamic Biosensors GmbH & Inc.  
BU-M-150-1 v2.1



## Product Description

Order Number: **BU-M-150-1**

*Table 1. Contents and Storage Information*

Material	Composition	Amount	Storage
<b>1x Buffer M pH 6.5</b>	50 mM 2-(N-Morpholino)ethanesulfonic acid, 150 mM NaCl; 0.2 µm sterile filtered	50 mL	2-8°C

For research use only.

This product has a limited shelf life, please see expiry date on label.

## Contact

**Dynamic Biosensors GmbH**  
Perchtinger Str. 8/10  
81379 Munich  
Germany

**Dynamic Biosensors, Inc.**  
300 Trade Center, Suite 1400  
Woburn, MA 01801  
USA

**Order Information** [order@dynamic-biosensors.com](mailto:order@dynamic-biosensors.com)  
**Technical Support** [support@dynamic-biosensors.com](mailto:support@dynamic-biosensors.com)

[www.dynamic-biosensors.com](http://www.dynamic-biosensors.com)

Instruments and chips are engineered and manufactured in Germany.  
©2024 Dynamic Biosensors GmbH | Dynamic Biosensors, Inc. All rights reserved.