# dynamic BIOSENSORS

# Safety Data Sheet

# according to UK REACH Regulation

# 10x HE140 Buffer, pH 7.4

Revision date: 07.12.2022

Product code: BU-HE-140-10

Page 1 of 9

# SECTION 1: Identification of the substance/mixture and of the company/undertaking

# 1.1. Product identifier

10x HE140 Buffer, pH 7.4

# 1.2. Relevant identified uses of the substance or mixture and uses advised against

# Use of the substance/mixture

Use as laboratory reagent The product is intended for research, analysis and scientific education.

Uses advised against

Any non-intended use.

# 1.3. Details of the supplier of the safety data sheet

Company name: Street: Place:	Dynamic Biosensors GmbH Perchtinger Str. 8/10 D-81379 München	
Telephone: Responsible Department:	+49-89-8974544-0 Dr. Gans-Eichler Chemieberatung GmbH Otto-Hahn-Str. 36 D-48161 Münster	e-mail: info@tge-consult.de Tel.: +49(0)2534 6441185 www.tge-consult.de
1.4. Emergency telephone	Emergency CONTACT (24-Ho	our-Number): GBK GmbH +49 (0)6132-84463

# 1.4. Emergency telephone

number:

# **Further Information**

Safety Data Sheet according to UK-REACH Regulation

# **SECTION 2: Hazards identification**

# 2.1. Classification of the substance or mixture

# **GB CLP Regulation**

This mixture is not classified as hazardous in accordance with GB CLP Regulation.

# 2.2. Label elements

# Additional advice on labelling

Labelling according to GHS (UK CLP) regulation .: none

# 2.3. Other hazards

The substances in the mixture (>0,1%) do not meet the PBT/vPvB criteria according to UK REACH This product does not contain a substance (> 0,1 %) that has endocrine disrupting properties with respect to non-target organisms as no components meets the criteria. No risks worthy of mention. Please observe the information on the safety data sheet at all times.

# **SECTION 3: Composition/information on ingredients**

# 3.2. Mixtures

# Chemical characterization

The product does not contain dangerous substances according to UK REACH, Annex II, Part A , 3.1/3.2. that must be mentioned in Chapter 3.

# **Further Information**

Product does not contain listed SVHC substances > 0,1 % according to UK REACH.

# **SECTION 4: First aid measures**

# 10x HE140 Buffer, pH 7.4

Revision date: 07.12.2022

Product code: BU-HE-140-10

Page 2 of 9

# 4.1. Description of first aid measures

# **General information**

In case of accident or unwellness, seek medical advice immediately (show directions for use or safety data sheet if possible).

# After inhalation

In case of accident by inhalation: remove casualty to fresh air and keep at rest. In case of respiratory tract irritation, consult a physician.

# After contact with skin

Gently wash with plenty of soap and water. In case of skin irritation, seek medical treatment.

#### After contact with eyes

Rinse cautiously with water for several minutes. In case of troubles or persistent symptoms, consult an ophthalmologist.

# After ingestion

Rinse mouth thoroughly with water. Let water be drunken in little sips (dilution effect). Do NOT induce vomiting. In all cases of doubt, or when symptoms persist, seek medical advice.

# 4.2. Most important symptoms and effects, both acute and delayed

No information available.

# 4.3. Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

# **SECTION 5: Firefighting measures**

# 5.1. Extinguishing media

#### Suitable extinguishing media

Carbon dioxide (CO2). Dry extinguishing powder. alcohol resistant foam. Atomized water.

# Unsuitable extinguishing media

High power water jet.

#### 5.2. Special hazards arising from the substance or mixture

Can be released in case of fire: Carbon monoxide Carbon dioxide (CO2). Hydrogen chloride (HCI), Nitrogen oxides (NOx), Sulphur oxides

# 5.3. Advice for firefighters

In case of fire: Wear self-contained breathing apparatus.

# Additional information

Collect contaminated fire extinguishing water separately. Do not allow entering drains or surface water. Co-ordinate fire-fighting measures to the fire surroundings.

# **SECTION 6: Accidental release measures**

# 6.1. Personal precautions, protective equipment and emergency procedures

# General advice

Safe handling: see section 7

#### For non-emergency personnel

Wear personal protection equipment (refer to section 8).

#### For emergency responders

No special measures are necessary.

# 6.2. Environmental precautions

Discharge into the environment must be avoided.

# 6.3. Methods and material for containment and cleaning up

#### For containment

Absorb with liquid-binding material (e.g. sand, diatomaceous earth, acid- or universal binding agents).

# dynamic BIOSENSORS

# Safety Data Sheet

# according to UK REACH Regulation

# 10x HE140 Buffer, pH 7.4

Revision date: 07.12.2022

Product code: BU-HE-140-10

Page 3 of 9

Treat the recovered material as prescribed in the section on waste disposal.

# For cleaning up

Clean contaminated objects and areas thoroughly observing environmental regulations.

# 6.4. Reference to other sections

Safe handling: see section 7 Personal protection equipment: see section 8 Disposal: see section 13

# **SECTION 7: Handling and storage**

# 7.1. Precautions for safe handling

#### Advice on safe handling

Wear suitable protective clothing. (See section 8.)

## Advice on protection against fire and explosion

Usual measures for fire prevention.

# Advice on general occupational hygiene

Always close containers tightly after the removal of product. Do not eat, drink, smoke or sneeze at the workplace. Wash hands before breaks and after work.

# Further information on handling

General protection and hygiene measures: See section 8.

#### 7.2. Conditions for safe storage, including any incompatibilities

#### Requirements for storage rooms and vessels

Keep container tightly closed in a cool, well-ventilated place.

#### Hints on joint storage

Do not store together with: Explosives. Oxidizing solids. Oxidizing liquids. Radioactive substances.. Infectious substances.

# Further information on storage conditions

Keep the packing dry and well sealed to prevent contamination and absorbtion of humidity. Recommended storage temperature: 2-8°C Protect against: frost. UV-radiation/sunlight. heat. Humidity

# 7.3. Specific end use(s)

See section 1.

#### **SECTION 8: Exposure controls/personal protection**

# 8.1. Control parameters

Additional advice on limit values

To date, no national critical limit values exist.

# 8.2. Exposure controls

# Appropriate engineering controls

Technical measures and the application of suitable work processes have priority over personal protection equipment.

Provide adequate ventilation.

# Individual protection measures, such as personal protective equipment

# Eye/face protection

Wear safety glasses; chemical goggles (if splashing is possible). BS/EN 166

# Hand protection

In case of prolonged or frequently repeated skin contact: Wear suitable gloves. Suitable material:

# 10x HE140 Buffer, pH 7.4

Revision date: 07.12.2022

Product code: BU-HE-140-10

Page 4 of 9

FKM (fluororubber). - Thickness of glove material: 0,4 mm Breakthrough time >= 8 h Butyl rubber. - Thickness of glove material: 0,5 mm Breakthrough time >= 8 h CR (polychloroprenes, Chloroprene rubber). - Thickness of glove material: 0,5 mm Breakthrough time >= 8 h NBR (Nitrile rubber). - Thickness of glove material: 0,35 mm Breakthrough time >= 8 h PVC (Polyvinyl chloride). - Thickness of glove material: 0,5 mm Breakthrough time >= 8 h The selected protective gloves have to satisfy the specifications of EU Directive EC/2016/425 and the standard EN ISO 374 derived from it. Before using check leak tightness / impermeability. In the case of wanting to use the gloves again, clean them before taking off and air them well.

# Skin protection

Suitable protective clothing: Lab apron.

Minimum standard for preventive measures while handling with working materials are specified in the TRGS 500 (D).

# **Respiratory protection**

With correct and proper use, and under normal conditions, breathing protection is not required.

#### Thermal hazards

No special precautionary measures are necessary.

# **Environmental exposure controls**

No special precautionary measures are necessary.

#### **SECTION 9: Physical and chemical properties**

# 9.1. Information on basic physical and chemical properties

1. Information on basic physical and chemical properties					
Physical state:	liquid				
Colour:	colourless				
Odour:	characteristic				
Odour threshold:	not determined				
Changes in the physical state					
Melting point/freezing point:		not determined			
Boiling point or initial boiling point and		not determined			
boiling range:					
Sublimation point:		not determined			
Softening point:		not determined			
Pour point:		not determined			
Flash point:		not determined			
Flammability					
Solid/liquid:		not determined			
Gas:		not relevant			
Explosive properties					
none					
Lower explosion limits:		not determined			
Upper explosion limits:		not determined			
Auto-ignition temperature:		not determined			
Self-ignition temperature					
Solid:		not relevant			

# dynamic BIOSENSORS

# Safety Data Sheet

# according to UK REACH Regulation

10x	HE140 Buffer, pH 7.4	
Revision date: 07.12.2022 Pro	duct code: BU-HE-140-10	Page 5 of 9
Gas:	not relevant	
Decomposition temperature:	not determined	
pH-Value:	7,4	
Viscosity / dynamic:	not determined	
Viscosity / kinematic:	not determined	
Flow time:	not determined	
Water solubility:	not determined	
Solubility in other solvents not determined Dissolution rate:	not relevant	
Partition coefficient n-octanol/water:	SECTION 12: Ecological information	
Dispersion stability:	not relevant	
Vapour pressure:	not determined	
Density: Bulk density:	not determined not determined	
Relative vapour density:	not determined	
Particle characteristics:	not relevant	
9.2. Other information		
Information with regard to physical hazard classes Sustaining combustion: Oxidizing properties none	Not sustaining combustion	
Other safety characteristics		
Solvent separation test:	not determined	
Solvent content:	not determined	
Solid content:	not determined	
Evaporation rate:	not determined	
Further Information No information available.		

SECTION 10: Stability and reactivity

# 10.1. Reactivity

No information available.

# 10.2. Chemical stability

The product is chemically stable under recommended conditions of storage, use and temperature.

# 10.3. Possibility of hazardous reactions

No hazardous reaction when handled and stored according to provisions. Refer to chapter 10.5.

# 10.4. Conditions to avoid

Protect against: UV-radiation/sunlight. heat.

# 10.5. Incompatible materials

Materials to avoid: Oxidizing agents, strong. Reducing agents, strong.

# 10.6. Hazardous decomposition products

Does not decompose when used for intended uses.

# 10x HE140 Buffer, pH 7.4

Revision date: 07.12.2022

Product code: BU-HE-140-10

Page 6 of 9

# SECTION 11: Toxicological information

# 11.1. Information on hazard classes as defined in GB CLP Regulation

# Toxicocinetics, metabolism and distribution

No data available.

# Acute toxicity

Based on available data, the classification criteria are not met. No data available.

## Irritation and corrosivity

Based on available data, the classification criteria are not met. No data available .

# Sensitising effects

Based on available data, the classification criteria are not met. No data available .

# Carcinogenic/mutagenic/toxic effects for reproduction

Based on available data, the classification criteria are not met. No data available

# STOT-single exposure

Based on available data, the classification criteria are not met. No data available .

#### STOT-repeated exposure

Based on available data, the classification criteria are not met. No data available .

# Aspiration hazard

Based on available data, the classification criteria are not met. No data available .

# Specific effects in experiment on an animal

No data available.

# 11.2. Information on other hazards

# Endocrine disrupting properties

This product does not contain a substance (> 0,1 %) that has endocrine disrupting properties with respect to non-target organisms as no components meets the criteria.

# Other information

No data available.

# **SECTION 12: Ecological information**

# 12.1. Toxicity

No data available

12.2. Persistence and degradability

No data available

# 12.3. Bioaccumulative potential

No data available

# 12.4. Mobility in soil

No data available.

# 12.5. Results of PBT and vPvB assessment

The substances in the mixture do not meet the PBT/vPvB criteria according to UK REACH. The aforementioned statement applies to substances contained in the product with a minimum content of 0.1

%.

# 10x HE140 Buffer, pH 7.4

Revision date: 07.12.2022

Product code: BU-HE-140-10

Page 7 of 9

# 12.6. Endocrine disrupting properties

This product does not contain a substance that has endocrine disrupting properties with respect to non-target organisms as no components meets the criteria.

The aforementioned statement applies to substances contained in the product with a minimum content of 0.1 %.

# 12.7. Other adverse effects

No data available.

# Further information

Do not allow to enter into surface water or drains.

# **SECTION 13: Disposal considerations**

# 13.1. Waste treatment methods

# Disposal recommendations

Observe in addition any national regulations! Consult the local waste disposal expert about waste disposal. Non-contaminated packages may be recycled.

According to (EWC) European Waste Catalogue, allocation of waste identity numbers/waste descriptions must be carried out in a specific way for every industry and process.

Control report for waste code/ waste marking according to (EWC) European Waste Catalogue:

# List of Wastes Code - residues/unused products

160509 WASTES NOT OTHERWISE SPECIFIED IN THE LIST; gases in pressure containers and discarded chemicals; discarded chemicals other than those mentioned in 16 05 06, 16 05 07 or 16 05 08

#### List of Wastes Code - used product

160509 WASTES NOT OTHERWISE SPECIFIED IN THE LIST; gases in pressure containers and discarded chemicals; discarded chemicals other than those mentioned in 16 05 06, 16 05 07 or 16 05 08

# List of Wastes Code - contaminated packaging

150106 WASTE PACKAGING; ABSORBENTS, WIPING CLOTHS, FILTER MATERIALS AND PROTECTIVE CLOTHING NOT OTHERWISE SPECIFIED; packaging (including separately collected municipal packaging waste); mixed packaging

# Contaminated packaging

Handle contaminated packages in the same way as the substance itself.

# **SECTION 14: Transport information**

# Land transport (ADR/RID)

14.1. UN number or ID number:	No dangerous good in sense of these transport regulations.
14.2. UN proper shipping name:	No dangerous good in sense of these transport regulations.
14.3. Transport hazard class(es):	No dangerous good in sense of these transport regulations.
14.4. Packing group:	No dangerous good in sense of these transport regulations.
Inland waterways transport (ADN)	
14.1. UN number or ID number:	No dangerous good in sense of these transport regulations.
14.2. UN proper shipping name:	No dangerous good in sense of these transport regulations.
14.3. Transport hazard class(es):	No dangerous good in sense of these transport regulations.
14.4. Packing group:	No dangerous good in sense of these transport regulations.
Marine transport (IMDG)	
14.1. UN number or ID number:	No dangerous good in sense of these transport regulations.
14.2. UN proper shipping name:	No dangerous good in sense of these transport regulations.
14.3. Transport hazard class(es):	No dangerous good in sense of these transport regulations.

Ir

N

	· ·	
Revision date: 07.12.2022	<b>10x HE140 Buffer, pH 7.4</b> Product code: BU-HE-140-10	Dage 9 of 0
Revision date. 07.12.2022		Page 8 of 9
Air transport (ICAO-TI/IATA-DGR)		
14.1. UN number or ID number:	No dangerous good in sense of these transport regulations.	
14.2. UN proper shipping name:	No dangerous good in sense of these transport regulations.	
14.3. Transport hazard class(es):	No dangerous good in sense of these transport regulations.	
14.5. Environmental hazards		
ENVIRONMENTALLY HAZARDOUS:	No	
14.6. Special precautions for user		
Refer to section 6 - 8		
14.7. Maritime transport in bulk according t	o IMO instruments	
not relevant		
SECTION 15: Regulatory information		
15.1. Safety, health and environmental regu	lations/legislation specific for the substance or mixture	
EU regulatory information		
Restrictions on use (REACH, annex XVII)		
Entry 75		
2010/75/EU (VOC):	No information available.	
2004/42/EC (VOC):	No information available.	
Information according to 2012/18/EU (SEVESO III):	Not subject to 2012/18/EU (SEVESO III)	
Additional information		
Safety Data Sheet according to UK-RI The mixture is classified as not hazard UK REACH Appendix XVII, No (mixtur	lous according to regulation (EC) No 1272/2008 [CLP].	
National regulatory information		
Water hazard class (D):	1 - slightly hazardous to water	
15.2. Chemical safety assessment		
For the following substances of this m	ixture a chemical safety assessment has been carried out:	
SECTION 16: Other information		
Changes		
Rev. 1.0; 08.02.2017, Initial release		
Rev. 2.0; 02.07.2021, Documentation Rev. 3.0; 07.12.2022; Changes in cha	of changes: chapter: 3, 6, 7, 8, 10, 11, 12, 15 ,16	
Abbreviations and acronyms		
-	t des marchandises dangereuses par Route (European Agreement	
concerning the International Carriage		
AGW: Arbeitsplatzgrenzwert		
CAS: Chemical Abstracts Service		
CLP: Classification, Labelling and Pac DNEL: Derived No Effect Level	kaging of substances and mixtures	
d: day(s)		
EINECS: European INventory of Exist	ing Commercial chemical Substances	
ELINCS: European List of Notified Ch	-	
ECHA: European Chemicals Agency		
EWC: European Waste Catalogue		
IARC: INTERNATIONAL AGENCY FC		

IMDG: International Maritime Code for Dangerous Goods

Page 9 of 9



# Safety Data Sheet

according to UK REACH Regulation

# **10x HE140 Buffer, pH 7.4** Product code: BU-HE-140-10

Revision date: 07.12.2022

IATA: International Air Transport Association IATA-DGR: Dangerous Goods Regulations by the "International Air Transport Association" (IATA) ICAO: International Civil Aviation Organization ICAO-TI: Technical Instructions by the "International Civil Aviation Organization" (ICAO) GHS: Globally Harmonized System of Classification and Labelling of Chemicals GefStoffV: Gefahrstoffverordnung (Ordinance on Hazardous Substances, Germany) h: hour LOAEL: Lowest observed adverse effect level LOAEC: Lowest observed adverse effect concentration LC50: Lethal concentration, 50 percent LD50: Lethal dose, 50 percent NOAEL: No observed adverse effect level NOAEC: No observed adverse effect concentration NLP: No-Longer Polymers N/A: not applicable OECD: Organisation for Economic Co-operation and Development PNEC: predicted no effect concentration PBT: Persistent bioaccumulative toxic RID: Reglement international concernant le transport des marchandises dangereuses par chemin de fer (Regulations Concerning the International Transport of Dangerous Goods by Rail) REACH: Registration, Evaluation, Authorisation of Chemicals SVHC: substance of very high concern TRGS: Technische Regeln für Gefahrstoffe UN: United Nations VOC: Volatile Organic Compounds

# **Further Information**

Classification according to GHS [UK CLP] - Classification procedure: Health hazards: Calculation method. Environmental hazards: Calculation method.

Physical hazards: On basis of test data and / or calculated and / or estimated.

The above information describes exclusively the safety requirements of the product and is based on our present-day knowledge. The information is intended to give you advice about the safe handling of the product named in this safety data sheet, for storage, processing, transport and disposal. The information cannot be transferred to other products. In the case of mixing the product with other products or in the case of processing, the information on this safety data sheet is not necessarily valid for the new made-up material.

(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)