

1x MES Buffer pH 6.5

This buffer is intended for use as a coupling buffer for nanolever conjugations of proteins with a low pI.

Product Description

Order Number **BU-M-150-1**

Material	Amount	Storage	Comment
1x MES Buffer pH 6.5 (50 mM 2-(N-Morpholino)ethanesulfonic acid, 150 mM NaCl)	12 mL	2 - 8°C	

For *in vitro* use only.

This product has a limited stability, please see expiry date on label.

Contact

Dynamic Biosensors GmbH

Perchtinger Str. 8/10
81379 Munich
Germany

Phone: +49 89 89 74 544 0

Dynamic Biosensors Inc.

300 Trade Center, Suite 1400
Woburn, MA 01801
USA

Phone: +1 781 404 6126

Order Information


Phone: +49 89 89 74 544 0

Email: order@dynamic-biosensors.com

Technical Support

Phone: +49 89 89 74 544 66

Email: support@dynamic-biosensors.com



switchSENSE® is a proprietary measurement technology by Dynamic Biosensors GmbH. Instruments and biochips are engineered and manufactured in Germany.

©2023 Dynamic Biosensors GmbH | Dynamic Biosensors Inc. All rights reserved.