

Safety Data Sheet

according to 29 CFR 1910.1200(g)

10 x PE40 pH 7.4

Revision date: 12/08/2022

Product code: BU-PE-40-10.

Page 1 of 9

1. Identification

Product identifier

10 x PE40 pH 7.4

Recommended use of the chemical and restrictions on use

Use of the substance/mixture

Use as laboratory reagent
 The product is intended for research, analysis and scientific education.

Uses advised against

Any non-intended use.

Details of the supplier of the safety data sheet

Company name:	Dynamic Biosensors Inc.	
Street:	300 Trade Center, Suite 1400	
Place:	USA-01801 Woburn, MA	
Telephone:	+1 781 404 6126	
Responsible Department:	Dr. Gans-Eichler	e-mail: info@tge-consult.de
	Chemieberatung GmbH	Tel.: +49(0)2534 6441185
	Otto-Hahn-Str. 36	www.tge-consult.de
	D-48161 Münster	

Emergency phone number: CONTACT (24-Hour-Number): GBK GmbH 01149-6132-84463

2. Hazard(s) identification

Classification of the chemical

29 CFR Part 1910.1200

Specific target organ toxicity repeated or prolonged exposure: STOT RE 2

Label elements

29 CFR Part 1910.1200

Signal word: Warning

Pictograms:



Hazard statements

May cause damage to organs through prolonged or repeated exposure

Precautionary statements

Do not breathe dust/fume/gas/mist/vapors/spray.
 Get medical advice/attention if you feel unwell.
 Dispose of contents/container to local/regional/national/international regulations.

Hazards not otherwise classified

The components in this formulation (>0,1%) do not meet the criteria for classification as PBT or vPvB.

3. Composition/information on ingredients

Mixtures

Hazardous components

CAS No	Components	Quantity
139-33-3	disodium dihydrogen ethylenediaminetetraacetate	2.5 %

Safety Data Sheet

according to 29 CFR 1910.1200(g)

10 x PE40 pH 7.4

Revision date: 12/08/2022

Product code: BU-PE-40-10.

Page 2 of 9

4. First-aid measures

Description of first aid measures

General information

In case of accident or unwellness, seek medical advice immediately (show directions for use or safety data sheet if possible).

After inhalation

In case of accident by inhalation: remove casualty to fresh air and keep at rest. In case of respiratory tract irritation, consult a physician.

After contact with skin

Gently wash with plenty of soap and water. In case of skin irritation, seek medical treatment.

After contact with eyes

Rinse cautiously with water for several minutes. In case of troubles or persistent symptoms, consult an ophthalmologist.

After ingestion

Rinse mouth thoroughly with water. Let water be drunk in little sips (dilution effect). Do NOT induce vomiting. In all cases of doubt, or when symptoms persist, seek medical advice.

Most important symptoms and effects, both acute and delayed

No information available.

Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

5. Fire-fighting measures

Extinguishing media

Suitable extinguishing media

Carbon dioxide (CO₂) Dry extinguishing powder. alcohol resistant foam. Atomized water.

Unsuitable extinguishing media

High power water jet.

Specific hazards arising from the chemical

Can be released in case of fire: Carbon monoxide. Carbon dioxide (CO₂). Phosphorus oxides. Nitrogen oxides (NO_x).

Special protective equipment and precautions for fire-fighters

In case of fire: Wear self-contained breathing apparatus.

Additional information

Collect contaminated fire extinguishing water separately. Do not allow entering drains or surface water.

Co-ordinate fire-fighting measures to the fire surroundings.

6. Accidental release measures

Personal precautions, protective equipment and emergency procedures

General advice

See protective measures under point 7 and 8.

For non-emergency personnel

Wear personal protection equipment (refer to section 8).

For emergency responders

No special measures are necessary.

Environmental precautions

Discharge into the environment must be avoided.

Methods and material for containment and cleaning up

Safety Data Sheet

according to 29 CFR 1910.1200(g)

10 x PE40 pH 7.4

Revision date: 12/08/2022

Product code: BU-PE-40-10.

Page 3 of 9

For containment

Absorb with liquid-binding material (e.g. sand, diatomaceous earth, acid- or universal binding agents).
Treat the recovered material as prescribed in the section on waste disposal.

For cleaning up

Clean contaminated objects and areas thoroughly observing environmental regulations.

Reference to other sections

Safe handling: see section 7
Personal protection equipment: see section 8
Disposal: see section 13

7. Handling and storage

Precautions for safe handling

Advice on safe handling

Wear suitable protective clothing. (See section 8.)

Advice on protection against fire and explosion

Usual measures for fire prevention.

Advice on general occupational hygiene

Always close containers tightly after the removal of product. Do not eat, drink, smoke or sneeze at the workplace. Wash hands before breaks and after work.

Further information on handling

General protection and hygiene measures: See section 8.

Conditions for safe storage, including any incompatibilities

Requirements for storage rooms and vessels

Keep container tightly closed in a cool, well-ventilated place.

Hints on joint storage

Do not store together with: Explosives. Oxidizing solids. Oxidizing liquids. Radioactive substances. Infectious substances. Food and animal feedingstuff.

Further information on storage conditions

Keep the packing dry and well sealed to prevent contamination and absorption of humidity.
Recommended storage temperature: 20°C
Protect against: frost. UV-radiation/sunlight. heat. Humidity

8. Exposure controls/personal protection

Control parameters

Additional advice on limit values

The following constituents are the only constituents of the product which have a PEL, TLV or other recommended exposure limit. At this time, the other constituents have no known exposure limits.

Exposure controls

Appropriate engineering controls

Technical measures and the application of suitable work processes have priority over personal protection equipment.
Provide adequate ventilation.

Individual protection measures, such as personal protective equipment

Eye/face protection

Wear safety glasses; chemical goggles (if splashing is possible). Standards: EN 166 or 29 CFR 1910.133

Hand protection

Wear suitable gloves.
Suitable material:
FKM (fluororubber). - Thickness of the glove material 0,4 mm

Safety Data Sheet

according to 29 CFR 1910.1200(g)

10 x PE40 pH 7.4

Revision date: 12/08/2022

Product code: BU-PE-40-10.

Page 4 of 9

Breakthrough time \geq 8 h

Butyl rubber. - Thickness of the glove material 0,5 mm

Breakthrough time \geq 8 h

CR (polychloroprenes, Chloroprene rubber). - Thickness of the glove material 0,5 mm

Breakthrough time \geq 8 h

NBR (Nitrile rubber). - Thickness of the glove material 0,35 mm

Breakthrough time \geq 8 h

PVC (Polyvinyl chloride). - Thickness of the glove material 0,5 mm

Breakthrough time \geq 8 h

For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves.

The selected protective gloves should satisfy the specifications of standards like EN 374.

Before using check leak tightness / impermeability. In the case of wanting to use the gloves again, clean them before taking off and air them well.

Skin protection

Suitable protective clothing: Lab apron.

Respiratory protection

With correct and proper use, and under normal conditions, breathing protection is not required.

Respiratory protection necessary at:

Exceeding exposure limit values

Suitable respiratory protective equipment: half-mask with filter EN 149 or 29 CFR 1910.134 .

The filter class must be suitable for the maximum contaminant concentration (gas/vapour/aerosol/particulates) that may arise when handling the product. If the concentration is exceeded, self-contained breathing apparatus must be used.

Environmental exposure controls

No special precautionary measures are necessary.

9. Physical and chemical properties

Information on basic physical and chemical properties

Physical state:	liquid
Color:	colourless
Odor:	odourless

Changes in the physical state

Melting point/freezing point:	not determined
Boiling point or initial boiling point and boiling range:	not determined
Sublimation point:	not determined
Softening point:	not determined
Pour point:	not determined
Flash point:	not determined

Explosive properties

none

Lower explosion limits:	not determined
Upper explosion limits:	not determined
Auto-ignition temperature:	not determined

Self-ignition temperature

Gas:	not determined
------	----------------

Decomposition temperature:	not determined
----------------------------	----------------

pH-Value (at 20 °C):	7,4
----------------------	-----

Safety Data Sheet

according to 29 CFR 1910.1200(g)

10 x PE40 pH 7.4

Revision date: 12/08/2022

Product code: BU-PE-40-10.

Page 5 of 9

Viscosity / dynamic: not determined
 Viscosity / kinematic: not determined
 Flow time: not determined
 Water solubility: completely miscible

Solubility in other solvents
 not determined

Partition coefficient n-octanol/water: SECTION 12: Ecological information

Vapor pressure: 23 hPa
 (at 20 °C)

Density (at 20 °C): 1,0277 g/cm³

Relative vapour density: not determined

Other information

Information with regard to physical hazard classes

Sustaining combustion: Not sustaining combustion

Oxidizing properties
 none

Other safety characteristics

Solvent separation test: not determined

Solvent content: Water: 94,9%

Solid content: not determined

Evaporation rate: not determined

Further Information

10. Stability and reactivity

Reactivity

No information available.

Chemical stability

Stability: Stable

The product is chemically stable under recommended conditions of storage, use and temperature.

Possibility of hazardous reactions

Hazardous reactions: Will not occur

No information available.

Conditions to avoid

Protect against: UV-radiation/sunlight. heat.

Incompatible materials

Materials to avoid: Oxidising agent, strong. Reducing agents, strong.

Hazardous decomposition products

Does not decompose when used for intended uses.

Can be released in case of fire: Carbon monoxide. Carbon dioxide (CO₂). Phosphorus oxides. Nitrogen oxides (NO_x).

11. Toxicological information

Route(s) of Entry

Ingestion: May be harmful if swallowed. Inhalation: May be harmful if inhaled. Skin contact: May cause irritation. Eye contact: May cause irritation.

Information on toxicological effects

Safety Data Sheet

according to 29 CFR 1910.1200(g)

10 x PE40 pH 7.4

Revision date: 12/08/2022

Product code: BU-PE-40-10.

Page 6 of 9

Toxicokinetics, metabolism and distribution

No data available.

Acute toxicity

Based on available data, the classification criteria are not met.

CAS No	Components				
	Exposure route	Dose	Species	Source	Method
139-33-3	disodium dihydrogen ethylenediaminetetraacetate				
	oral	LD50 2800 mg/kg	Rat	ECHA Dossier	
	inhalation vapour	ATE 11 mg/l			
	inhalation dust/mist	ATE 1,5 mg/l			

Irritation and corrosivity

Based on available data, the classification criteria are not met.

Sensitizing effects

Based on available data, the classification criteria are not met.

Carcinogenic/mutagenic/toxic effects for reproduction

Based on available data, the classification criteria are not met.

Specific target organ toxicity (STOT) - single exposure

Based on available data, the classification criteria are not met.

Specific target organ toxicity (STOT) - repeated exposure

May cause damage to organs through prolonged or repeated exposure (disodium dihydrogen ethylenediaminetetraacetate)

Carcinogenicity (OSHA): No ingredient of this mixture is listed.

Carcinogenicity (IARC): No ingredient of this mixture is listed.

Carcinogenicity (NTP): No ingredient of this mixture is listed.

Aspiration hazard

Based on available data, the classification criteria are not met.

Specific effects in experiment on an animal

No data available.

Information on other hazards

Endocrine disrupting properties

No data available.

12. Ecological information

Ecotoxicity

The product has not been tested.

Persistence and degradability

The product has not been tested.

Bioaccumulative potential

No indication of bioaccumulation potential.

Mobility in soil

No data available.

Endocrine disrupting properties

This product does not contain a substance that has endocrine disrupting properties with respect to non-target organisms as no components meets the criteria.

The aforementioned statement applies to substances contained in the product with a minimum content of 0.1%.

Safety Data Sheet

according to 29 CFR 1910.1200(g)

10 x PE40 pH 7.4

Revision date: 12/08/2022

Product code: BU-PE-40-10.

Page 7 of 9

Other adverse effects

No data available.

Further information

Do not allow to enter into surface water or drains.

13. Disposal considerations

Waste treatment methods

Disposal recommendations

Observe in addition any national regulations! Consult the local waste disposal expert about waste disposal.
Non-contaminated packages may be recycled.

Contaminated packaging

Handle contaminated packages in the same way as the substance itself.

14. Transport information

US DOT 49 CFR 172.101

Proper shipping name:

Not a hazardous material with respect to these transport regulations. &&
Not controlled under DOT

Marine transport (IMDG)

UN number or ID number:

No dangerous good in sense of this transport regulation.

UN proper shipping name:

No dangerous good in sense of this transport regulation.

Transport hazard class(es):

No dangerous good in sense of this transport regulation.

Packing group:

No dangerous good in sense of this transport regulation.

Air transport (ICAO-TI/IATA-DGR)

UN number or ID number:

No dangerous good in sense of this transport regulation.

UN proper shipping name:

No dangerous good in sense of this transport regulation.

Transport hazard class(es):

No dangerous good in sense of this transport regulation.

Packing group:

No dangerous good in sense of this transport regulation.

Environmental hazards

ENVIRONMENTALLY HAZARDOUS: No

Special precautions for user

See section 8.

Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

not relevant

15. Regulatory information

U.S. Regulations

National Inventory TSCA

disodium dihydrogen ethylenediaminetetraacetate listed in the TSCA inventory 8 (b): (x) active ,not listed under TSCA 12(b)

National regulatory information

SARA Section 311/312 Hazards:

disodium dihydrogen ethylenediaminetetraacetate (139-33-3): Immediate (acute) health hazard, Delayed (chronic) health hazard

State Regulations

Safe Drinking Water and Toxic Enforcement Act of 1986 (Proposition 65, State of California)

This product can not expose you to chemicals known to the State of California to cause cancer, birth defects or

Safety Data Sheet

according to 29 CFR 1910.1200(g)

10 x PE40 pH 7.4

Revision date: 12/08/2022

Product code: BU-PE-40-10.

Page 8 of 9

other reproductive harm.

This preparation is hazardous in the sense of regulation 29 CFR Part 1910.1200.

16. Other information

Hazardous Materials Information Label (HMIS)

Health: *1
 Flammability: 0
 Physical Hazard: 0
 Personal Protection: -

NFPA Hazard Ratings

Health: 1
 Flammability: 0
 Reactivity: 0
 Unique Hazard: -



Changes

Revision date: 08.12.2022
 Revision No: 2,0
 Rev. 1.0; Initial release: 05.07.2021
 Rev. 2.0; 08.12.2022; Changes in chapter: 1,16

Abbreviations and acronyms

- ACGIH: American Conference of Governmental Industrial Hygienists
- ASTM: American Society for Testing and Materials.
- ATE: acute toxicity estimate
- BCF: Bio concentration factor
- ECHA: European Chemicals Agency
- CAS: Chemical Abstracts Service
- CFR: Code of Federal Regulations
- DOT: Department of Transportation
- d: days
- EC50: Half maximal effective concentration
- EN: European Norm
- EPA: Environmental Protection Agency
- GHS: Globally Harmonized System of Classification and Labelling of Chemicals
- h: hours
- HMIS: Hazardous Materials Identification System
- IARC: INTERNATIONAL AGENCY FOR RESEARCH ON CANCER
- IBC: Intermediate Bulk Container
- IMDG: International Maritime Code for Dangerous Goods
- IATA: International Air Transport Association
- IATA-DGR: Dangerous Goods Regulations by the "International Air Transport Association" (IATA)
- ICAO: International Civil Aviation Organization
- ICAO-TI: Technical Instructions by the "International Civil Aviation Organization" (ICAO)
- GHS: Globally Harmonized System of Classification and Labelling of Chemicals
- LOAEL: Lowest observed adverse effect level
- LOAEC: Lowest observed adverse effect concentration
- LC50: Lethal concentration, 50 percent
- LD50: Lethal dose, 50 percent
- MARPOL: marine pollution
- NOAEL: No observed adverse effect level
- NOAEC: No observed adverse effect concentration
- NTP: National Toxicology Program
- N/A: not applicable

Safety Data Sheet

according to 29 CFR 1910.1200(g)

10 x PE40 pH 7.4

Revision date: 12/08/2022

Product code: BU-PE-40-10.

Page 9 of 9

NFPA: National Fire Protection Association
UN: United Nations
OECD: Organisation for Economic Co-operation and Development
OSHA: Occupational Safety and Health Administration
PBT: Persistent bioaccumulative toxic
RTECS: Registry of Toxic Effects of Chemical Substances
REACH: Registration, Evaluation, Authorisation and Restriction of Chemicals
SARA: Superfund Amendments and Reauthorization Act
STEL: short-term exposure limits
TSCA: Toxic Substances Control Act
TWA: time weighted average
VOC: Volatile Organic Compounds

Other data

Classification according 29 CFR Part 1910.1200: - Classification procedure:
Health hazards: Calculation method.
Environmental hazards: Calculation method.
Physical hazards: On basis of test data and / or calculated and / or estimated.

The above information describes exclusively the safety requirements of the product and is based on our present-day knowledge. The information is intended to give you advice about the safe handling of the product named in this safety data sheet, for storage, processing, transport and disposal. The information cannot be transferred to other products. In the case of mixing the product with other products or in the case of processing, the information on this safety data sheet is not necessarily valid for the new made-up material.

(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)