

according to UK REACH Regulation

**TCEP** 

Revision date: 05.12.2022 Product code: RK-PA-1 Page 1 of 11

## SECTION 1: Identification of the substance/mixture and of the company/undertaking

#### 1.1. Product identifier

**TCEP** 

Substance name: Tris(2-carboxyethyl)phosphine hydrochloride

CAS No: 51805-45-9 EC No: 629-759-1

#### 1.2. Relevant identified uses of the substance or mixture and uses advised against

## Use of the substance/mixture

Use as laboratory reagent. Scientific research and development.

#### Uses advised against

Any non-intended use.

## 1.3. Details of the supplier of the safety data sheet

Company name: Dynamic Biosensors GmbH Street: Perchtinger Str. 8/10 Place: D-81379 München Telephone: +49-89-8974544-0

Responsible Department:

Dr. Gans-Eichler

Chemieberatung GmbH

Otto-Hahn-Str. 36

Dr. Gans-Eichler

e-mail: info@tge-consult.de

D-48161 Münster

1.4. Emergency telephone

Emergency CONTACT (24-Hour-Number): GBK GmbH +49 (0)6132-84463

number:

#### **Further Information**

Safety Data Sheet according to UK-REACH Regulation

#### **SECTION 2: Hazards identification**

## 2.1. Classification of the substance or mixture

#### **GB CLP Regulation**

Skin Corr. 1B; H314

Full text of hazard statements: see SECTION 16.

#### 2.2. Label elements

#### **GB CLP Regulation**

Signal word: Danger

Pictograms:



# **Hazard statements**

H314 Causes severe skin burns and eye damage.

## **Precautionary statements**

P264 Wash hands thoroughly after handling.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P301+P330+P331 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.

P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with

water or shower.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if



#### according to UK REACH Regulation

**TCEP** 

Revision date: 05.12.2022 Product code: RK-PA-1 Page 2 of 11

present and easy to do. Continue rinsing.

P501 Dispose of contents/container in accordance with local/regional/national/international

regulations.

#### 2.3. Other hazards

The substances in the mixture (>0,1%) do not meet the PBT/vPvB criteria according to UK REACH This product does not contain a substance (> 0,1 %) that has endocrine disrupting properties with respect to non-target organisms as no components meets the criteria.

## **SECTION 3: Composition/information on ingredients**

#### 3.1. Substances

Sum formula: C9H15O6P · HCI
Molecular weight: 286,65 g/mol

#### **Hazardous components**

CAS No	Chemical name			
	EC No	Index No	REACH No	
	GHS Classification			
51805-45-9	Tris(2-carboxyethyl)phosphine hydrochloride			99 - < 100 %
	629-759-1			
	Skin Corr. 1B; H314			

Full text of H and EUH statements: see section 16.

## Specific Conc. Limits, M-factors and ATE

CAS No	EC No	Chemical name	Quantity
	Specific Conc. Limits, M-factors and ATE		
51805-45-9	629-759-1	Tris(2-carboxyethyl)phosphine hydrochloride	99 - < 100 %
	dermal: LD50 = 3000 mg/kg; oral: LD50 = 3500 mg/kg		

#### **Further Information**

Product does not contain listed SVHC substances > 0,1 % according to UK REACH.

# **SECTION 4: First aid measures**

## 4.1. Description of first aid measures

## **General information**

In case of accident or unwellness, seek medical advice immediately (show directions for use or safety data sheet if possible).

#### After inhalation

In case of accident by inhalation: remove casualty to fresh air and keep at rest. In case of respiratory tract irritation, consult a physician. In the case of lung irritation: Primary treatment using corticoide spray, eg. Auxiloson spray, Pulmicort-dosage-spray. (Auxiloson and Pulmicort are registered trademarks).

## After contact with skin


After contact with skin, wash immediately with plenty of water and soap. Take off immediately all contaminated clothing. In case of skin irritation, seek medical treatment.

#### After contact with eyes

In case of contact with eyes, rinse immediately with plenty of flowing water for 10 to 15 minutes holding eyelids apart. Subsequently consult an ophthalmologist.

# After ingestion

Do NOT induce vomiting. Rinse mouth thoroughly with water. Let water be drunken in little sips (dilution effect). Observe risk of aspiration if vomiting occurs. Never give anything by mouth to an unconscious person or a person with cramps. When in doubt or if symptoms are observed, get medical advice.





according to UK REACH Regulation

#### **TCEP**

Revision date: 05.12.2022 Product code: RK-PA-1 Page 3 of 11

#### 4.2. Most important symptoms and effects, both acute and delayed

If swallowed danger of perforation of the esophagus and the stomach (strong corrosive effects).

## 4.3. Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

## **SECTION 5: Firefighting measures**

## 5.1. Extinguishing media

#### Suitable extinguishing media

Sand. Foam. Carbon dioxide (CO2). Extinguishing powder. In case of major fire and large quantities: Water spray jet. Water mist.

## Unsuitable extinguishing media

High power water jet.

# 5.2. Special hazards arising from the substance or mixture

Can be released in case of fire: Carbon monoxide Carbon dioxide (CO2)

#### 5.3. Advice for firefighters

In case of fire and/or explosion do not breathe fumes. In case of fire: Wear self-contained breathing apparatus.

#### Additional information

Collect contaminated fire extinguishing water separately. Do not allow entering drains or surface water.

Co-ordinate fire-fighting measures to the fire surroundings.

#### **SECTION 6: Accidental release measures**

## 6.1. Personal precautions, protective equipment and emergency procedures

## General advice

Do not breathe vapour/aerosol. Avoid contact with skin, eyes and clothes.

## For non-emergency personnel

Wear personal protection equipment (refer to section 8).

# For emergency responders

No special measures are necessary.

#### 6.2. Environmental precautions

Do not allow to enter into surface water or drains. Prevent spread over a wide area (e.g. by containment or oil barriers). Do not allow to enter into soil/subsoil.

#### 6.3. Methods and material for containment and cleaning up

#### For containment

Absorb with liquid-binding material (e.g. sand, diatomaceous earth, acid- or universal binding agents).

Treat the recovered material as prescribed in the section on waste disposal.

#### For cleaning up

Clean contaminated objects and areas thoroughly observing environmental regulations.

## 6.4. Reference to other sections

Safe handling: see section 7

Personal protection equipment: see section 8

Disposal: see section 13

## **SECTION 7: Handling and storage**

#### 7.1. Precautions for safe handling

#### Advice on safe handling

Wear suitable protective clothing. (See section 8.)

Conditions to avoid: aerosol or mist formation

Avoid contact with skin, eyes and clothes.



according to UK REACH Regulation

#### **TCEP**

Revision date: 05.12.2022 Product code: RK-PA-1 Page 4 of 11

#### Advice on protection against fire and explosion

Usual measures for fire prevention.

#### Advice on general occupational hygiene

When using do not eat, drink or smoke.

#### Further information on handling

General protection and hygiene measures: See section 8.

## 7.2. Conditions for safe storage, including any incompatibilities

#### Requirements for storage rooms and vessels

Keep container tightly closed in a cool, well-ventilated place. Only use containers specifically approved for the substance/product.

Make sure spills can be contained (e.g. sump pallets or kerbed areas).

#### Hints on joint storage

Do not store together with: Explosives. Oxidizing solids. Oxidizing liquids. Organic peroxides. Self-reactive substances and mixtures. Radioactive substances.. Infectious substances.

#### Further information on storage conditions

Recommended storage temperature: -20°C

Protect against: frost. UV-radiation/sunlight. heat. Humidity

#### 7.3. Specific end use(s)

See section 1.

## **SECTION 8: Exposure controls/personal protection**

#### 8.1. Control parameters

# Additional advice on limit values

To date, no national critical limit values exist.

## 8.2. Exposure controls







#### Appropriate engineering controls

Technical measures and the application of suitable work processes have priority over personal protection equipment.

Provide adequate ventilation.

#### Individual protection measures, such as personal protective equipment

# Eye/face protection

Wear eye/face protection. BS/EN 166

## **Hand protection**

Wear suitable gloves.

Suitable material:

FKM (fluororubber). - Thickness of glove material: 0,4 mm

Breakthrough time >= 8 h

Butyl rubber. - Thickness of glove material: 0,5 mm

Breakthrough time >= 8 h

CR (polychloroprenes, Chloroprene rubber). - Thickness of glove material: 0,5 mm

Breakthrough time >= 8 h

NBR (Nitrile rubber). - Thickness of glove material: 0,35 mm

Breakthrough time >= 8 h

PVC (Polyvinyl chloride). - Thickness of glove material: 0,5 mm

Breakthrough time >= 8 h



according to UK REACH Regulation

#### **TCEP**

Revision date: 05.12.2022 Product code: RK-PA-1 Page 5 of 11

For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves.

The selected protective gloves have to satisfy the specifications of EU Directive EC/2016/425 and the standard EN ISO 374 derived from it.

Before using check leak tightness / impermeability. In the case of wanting to use the gloves again, clean them before taking off and air them well.

#### Skin protection

Suitable protective clothing: Lab apron.

Minimum standard for preventive measures while handling with working materials are specified in the TRGS 500 (D).

#### Respiratory protection

With correct and proper use, and under normal conditions, breathing protection is not required.

Respiratory protection necessary at:

- -Exceeding exposure limit values
- -Insufficient ventilation and aerosol or mist formation

Suitable respiratory protective equipment: particulates filter device (DIN EN 143). Type: P1-3

The filter class must be suitable for the maximum contaminant concentration (gas/vapour/aerosol/particulates) that may arise when handling the product. If the concentration is exceeded, self-contained breathing apparatus must be used.

#### **Environmental exposure controls**

No information available.

# SECTION 9: Physical and chemical properties

#### 9.1. Information on basic physical and chemical properties

Physical state: solid Colour: white

Odour: characteristic
Odour threshold: not determined

#### Changes in the physical state

Melting point/freezing point:

No information available.

No information available.

No information available.

boiling range:

Sublimation point:No information available.Softening point:No information available.Pour point:No information available.Flash point:No information available.

Flammability

Solid/liquid: No information available.

Gas: not relevant

**Explosive properties** 

none

Lower explosion limits:

Upper explosion limits:

No information available.

No information available.

No information available.

Self-ignition temperature

Solid: not relevant Gas: not relevant Decomposition temperature: No information available. pH-Value: No information available.



according to UK REACH Regulation

**TCEP** 

Revision date: 05.12.2022 Product code: RK-PA-1 Page 6 of 11

Viscosity / dynamic:

Viscosity / kinematic:

No information available.

No information available.

No information available.

No information available.

Immiscible

Solubility in other solvents

No information available.

not relevant Dissolution rate: Partition coefficient n-octanol/water: No information available. Dispersion stability: not relevant Vapour pressure: No information available. Density: No information available. No information available. Bulk density: No information available. Relative vapour density: Particle characteristics: not relevant

9.2. Other information

Information with regard to physical hazard classes

Sustaining combustion: No data available

Oxidizing properties

none

Other safety characteristics

Solvent separation test:

Solvent content:

No information available.

Solid content:

No information available.

No information available.

Evaporation rate:

No information available.

Further Information

No information available.

# **SECTION 10: Stability and reactivity**

#### 10.1. Reactivity

No information available.

## 10.2. Chemical stability

The product is chemically stable under recommended conditions of storage, use and temperature.

## 10.3. Possibility of hazardous reactions

No hazardous reaction when handled and stored according to provisions.

Refer to chapter 10.5.

## 10.4. Conditions to avoid

Protect against: UV-radiation/sunlight. heat.

# 10.5. Incompatible materials

Materials to avoid: Oxidizing agents, strong. Reducing agents, strong.

## 10.6. Hazardous decomposition products

Does not decompose when used for intended uses.

## **SECTION 11: Toxicological information**

# 11.1. Information on hazard classes as defined in GB CLP Regulation



according to UK REACH Regulation

## **TCEP**

Revision date: 05.12.2022 Product code: RK-PA-1 Page 7 of 11

## Toxicocinetics, metabolism and distribution

No information available.

#### **Acute toxicity**

Based on available data, the classification criteria are not met.

CAS No	Chemical name						
	Exposure route	Dose	Species	Source	Method		
51805-45-9	Tris(2-carboxyethyl)phosphine hydrochloride						
	oral	LD50 3500 mg/kg	Rat				
	dermal	LD50 3000 mg/kg	Rat				

## Irritation and corrosivity

Causes severe skin burns and eye damage.

Causes serious eye damage.

#### Sensitising effects

Based on available data, the classification criteria are not met.

#### Carcinogenic/mutagenic/toxic effects for reproduction

Based on available data, the classification criteria are not met.

#### STOT-single exposure

Based on available data, the classification criteria are not met.

#### STOT-repeated exposure

Based on available data, the classification criteria are not met.

#### **Aspiration hazard**

Based on available data, the classification criteria are not met.

## 11.2. Information on other hazards

#### **Endocrine disrupting properties**

This product does not contain a substance (> 0,1 %) that has endocrine disrupting properties with respect to non-target organisms as no components meets the criteria.

# Other information

No data available.

## **SECTION 12: Ecological information**

#### 12.1. Toxicity

The product has not been tested.

## 12.2. Persistence and degradability

The product has not been tested.

#### 12.3. Bioaccumulative potential

No indication of bioaccumulation potential.

## 12.4. Mobility in soil

No information available.

## 12.5. Results of PBT and vPvB assessment

This substance does not meet the PBT/vPvB criteria of UK REACH.

The aforementioned statement applies to substances contained in the product with a minimum content of 0.1 %.

## 12.6. Endocrine disrupting properties

This substance does not have endocrine disrupting properties with respect to non-target organisms.

The aforementioned statement applies to substances contained in the product with a minimum content of 0.1 %



according to UK REACH Regulation

#### **TCEP**

Revision date: 05.12.2022 Product code: RK-PA-1 Page 8 of 11

#### 12.7. Other adverse effects

No information available.

## **SECTION 13: Disposal considerations**

#### 13.1. Waste treatment methods

#### **Disposal recommendations**

Dispose of waste according to applicable legislation. Consult the local waste disposal expert about waste disposal. Non-contaminated packages may be recycled. The allocation of waste identity numbers/waste descriptions must be carried out according to the EEC, specific to the industry and process. Control report for waste code/ waste marking according to (EWC) European Waste Catalogue:

#### List of Wastes Code - residues/unused products

160508 WASTES NOT OTHERWISE SPECIFIED IN THE LIST; gases in pressure containers and

discarded chemicals; discarded organic chemicals consisting of or containing hazardous

substances; hazardous waste

## List of Wastes Code - used product

160508 WASTES NOT OTHERWISE SPECIFIED IN THE LIST; gases in pressure containers and

discarded chemicals; discarded organic chemicals consisting of or containing hazardous

substances; hazardous waste

#### List of Wastes Code - contaminated packaging

150110 WASTE PACKAGING; ABSORBENTS, WIPING CLOTHS, FILTER MATERIALS AND

PROTECTIVE CLOTHING NOT OTHERWISE SPECIFIED; packaging (including separately collected municipal packaging waste); packaging containing residues of or contaminated by

hazardous substances; hazardous waste

## Contaminated packaging

Handle contaminated packages in the same way as the substance itself.

## **SECTION 14: Transport information**

#### Land transport (ADR/RID)

14.1. UN number or ID number: UN 3261

14.2. UN proper shipping name: CORROSIVE SOLID, ACIDIC, ORGANIC, N.O.S.

(Tris(2-carboxyethyl)phosphine hydrochloride)

14.3. Transport hazard class(es): 8

14.4. Packing group:

Hazard label: 8



Classification code: C4
Special Provisions: 274
Limited quantity: 1 kg
Excepted quantity: E2
Transport category: 2
Hazard No: 80
Tunnel restriction code: E

Inland waterways transport (ADN)

14.1. UN number or ID number: UN 3261

14.2. UN proper shipping name: CORROSIVE SOLID, ACIDIC, ORGANIC, N.O.S.

(Tris(2-carboxyethyl)phosphine hydrochloride)

14.3. Transport hazard class(es): 8

# **dynamic**BIOSENSORS

# **Safety Data Sheet**

according to UK REACH Regulation

**TCEP** 

Revision date: 05.12.2022 Product code: RK-PA-1 Page 9 of 11

14.4. Packing group:

Hazard label: 8



Ш

Classification code: C4
Special Provisions: 274
Limited quantity: 1 kg
Excepted quantity: E2

Marine transport (IMDG)

14.1. UN number or ID number: UN 3261

14.2. UN proper shipping name: CORROSIVE SOLID, ACIDIC, ORGANIC, N.O.S.

(Tris(2-carboxyethyl)phosphine hydrochloride)

14.3. Transport hazard class(es):814.4. Packing group:IIHazard label:8



Marine pollutant:

Special Provisions:

Limited quantity:

Excepted quantity:

EmS:

F-A, S-B

Segregation group:

NO

274

1 kg

E2

F-A, S-B

1 - acids

Air transport (ICAO-TI/IATA-DGR)

14.1. UN number or ID number: UN 3261

14.2. UN proper shipping name: CORROSIVE SOLID, ACIDIC, ORGANIC, N.O.S.

(Tris(2-carboxyethyl)phosphine hydrochloride)

14.3. Transport hazard class(es):814.4. Packing group:IIHazard label:8



Special Provisions:

Limited quantity Passenger:

Passenger LQ:

Excepted quantity:

A3 A803

5 kg

Y844

Excepted quantity:

E2

IATA-packing instructions - Passenger: 859
IATA-max. quantity - Passenger: 15 kg
IATA-packing instructions - Cargo: 863
IATA-max. quantity - Cargo: 50 kg

14.5. Environmental hazards

ENVIRONMENTALLY HAZARDOUS: No

14.6. Special precautions for user

Safe handling: see section 7

Personal protection equipment: see section 8



according to UK REACH Regulation

**TCEP** 

Revision date: 05.12.2022 Product code: RK-PA-1 Page 10 of 11

#### 14.7. Maritime transport in bulk according to IMO instruments

not relevant

# **SECTION 15: Regulatory information**

#### 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

#### **EU** regulatory information

Restrictions on use (REACH, annex XVII):

Entry 3

2010/75/EU (VOC):

No information available.

2004/42/EC (VOC):

No information available.

Information according to 2012/18/EU Not subject to 2012/18/EU (SEVESO III)

(SEVESO III):

#### **Additional information**

Safety Data Sheet according to UK-REACH Regulation

The mixture is classified as hazardous according to regulation (EC) No 1272/2008 [CLP].

UK REACH Appendix XVII, No (mixture): 3

#### **National regulatory information**

Employment restrictions: Observe restrictions to employment for juveniles according to the 'juvenile

work protection guideline' (94/33/EC).

Water hazard class (D): 3 - highly hazardous to water

#### 15.2. Chemical safety assessment

For the following substances of this mixture a chemical safety assessment has been carried out:

## **SECTION 16: Other information**

#### Changes

Rev. 1.00; Initial release: 05.12.2022

#### Abbreviations and acronyms

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement

concerning the International Carriage of Dangerous Goods by Road)

CAS: Chemical Abstracts Service

CLP: Classification, Labelling and Packaging of substances and mixtures

DNEL: Derived No Effect Level

d: day(s)

EINECS: European INventory of Existing Commercial chemical Substances

ELINCS: European List of Notified Chemical Substances

ECHA: European Chemicals Agency EWC: European Waste Catalogue

IARC: INTERNATIONAL AGENCY FOR RESEARCH ON CANCER

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

IATA-DGR: Dangerous Goods Regulations by the "International Air Transport Association" (IATA)

ICAO: International Civil Aviation Organization

ICAO-TI: Technical Instructions by the "International Civil Aviation Organization" (ICAO)

GHS: Globally Harmonized System of Classification and Labelling of Chemicals GefStoffV: Gefahrstoffverordnung (Ordinance on Hazardous Substances, Germany)


h: hour

LOAEL: Lowest observed adverse effect level

LOAEC: Lowest observed adverse effect concentration

LC50: Lethal concentration, 50 percent

LD50: Lethal dose, 50 percent





#### according to UK REACH Regulation

**TCEP** 

Revision date: 05.12.2022 Product code: RK-PA-1 Page 11 of 11

NOAEL: No observed adverse effect level

NOAEC: No observed adverse effect concentration

NLP: No-Longer Polymers

N/A: not applicable

OECD: Organisation for Economic Co-operation and Development

PNEC: predicted no effect concentration PBT: Persistent bioaccumulative toxic

RID: Règlement international concernant le transport des marchandises dangereuses par chemin de

fer (Regulations Concerning the International Transport of Dangerous Goods by Rail )

REACH: Registration, Evaluation, Authorisation of Chemicals

SVHC: substance of very high concern TRGS: Technische Regeln für Gefahrstoffe

**UN: United Nations** 

VOC: Volatile Organic Compounds

## Relevant H and EUH statements (number and full text)

H314 Causes severe skin burns and eye damage.

H318 Causes serious eye damage.

#### **Further Information**

Classification according to GHS [UK CLP] - Classification procedure:

Health hazards: Calculation method. Environmental hazards: Calculation method.

Physical hazards: On basis of test data and / or calculated and / or estimated.

The above information describes exclusively the safety requirements of the product and is based on our present-day knowledge. The information is intended to give you advice about the safe handling of the product named in this safety data sheet, for storage, processing, transport and disposal. The information cannot be transferred to other products. In the case of mixing the product with other products or in the case of processing, the information on this safety data sheet is not necessarily valid for the new made-up material.