

Urea solution

This solution is intended for use as an alternate regeneration solution for a **switch**SENSE® chip.

Product description

Order Number **SOL-REG-UR-12-1**

Material	Amount	Storage	Comment
Urea solution (7 M Urea; 0.2 µm sterile filtered)	12 mL	2-8°C	

For *in vitro* use only.

This product has a limited stability, please see expiry date on label.

Please be aware that this solution is not suitable for a complete regeneration of the chip as only proteins will be removed. The complementary nanolever will stay on the chip.

Preparation

This solution can be used without any further treatment.

Contact

Dynamic Biosensors GmbH

Perchtinger Str. 8/10
81379 Munich
Germany

Phone: +49 89 89 74 544 0

Dynamic Biosensors Inc.

300 Trade Center, Suite 1400
Woburn, MA 01801
USA

Phone: +1 781 404 6126

Order Information


Phone: +49 89 89 74 544 0

Email: order@dynamic-biosensors.com

Technical Support

Phone: +49 89 89 74 544 66

Email: support@dynamic-biosensors.com



switchSENSE® is a proprietary measurement technology by Dynamic Biosensors GmbH. Instruments and biochips are engineered and manufactured in Germany.

©2023 Dynamic Biosensors GmbH | Dynamic Biosensors Inc. All rights reserved.